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1/3

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1. Trapani JA: Granzymes: a family of lymphocyte granule serine proteases. *Genome Biol* 2: REVIEWS 3014, 2001.
2. Kam C-M et al., Granzymes (lymphocyte serine proteases): characterization with natural and synthetic substrates and inhibitors, *Biochimica et Biophysica Acta*, 2000, 1477, 307-323.
3. Estebanez-Perpina E, Fuentes-Prior P, Belorgey D, et al: Crystal structure of the caspase activator human granzyme B, a proteinase highly specific for an Asp-P1 residue. *Biol Chem* 381: 1203-1214, 2000.
4. Shi L, Keefe D, Durand E, Feng H, Zhang D and Lieberman J: Granzyme B binds to target cells mostly by charge and must be added at the same time as perforin to trigger apoptosis. *J Immunol* 174: 5456-5461, 2005.
5. Pinkoski MJ, Hobman M, Heibein JA, et al: Entry and trafficking of granzyme B in target cells during granzyme B perforin- mediated apoptosis. *Blood* 92: 1044-1054, 1998.
6. Motyka B, Korbutt G, Pinkoski MJ, et al: Mannose 6-phosphate/insulin-like growth factor II receptor is a death receptor for granzyme B during cytotoxic T cell-induced apoptosis. *Cell* 103: 491-500, 2000.
7. Jahrsdorfer B. et al, Granzyme B produced by human plasmacytoid dendritic cells suppresses T-cell expansion, *Blood*, 2010, 11, 115, 6.
8. Grifliths GM and Isaaz S., Granzymes A and B Are Targeted to the Lytic Granules of Lymphocytes by the Mannose-6-Phosphate Receptor, *The Journal of Cell Biology*, 1993, 120, 4, 885-896.
9. Thiery et al., Perforin activates clathrin- and dynamin-dependent endocytosis, which is required for plasma membrane repair and delivery of granzyme B for granzyme-mediated apoptosis, *Blood*, 2010, 115(8), 1582–1593.
10. Wagner et al., Expression of granzyme B in peripheral blood polymorpho-nuclear neutrophils (PMN), myeloid cell lines and in PMN derived from haematopoietic stem cells in vitro, *Molecular Immunology* 2008, 45, 1761–1766.

免疫泰克有限公司 IMMUNOTECH S.A.S. (a Beckman Coulter Company) , 130, avenue de Lattre de Tassigny, BP 177, 13276 Marseille cedex 9, France, 33-491 172 727